

ECIIS
EUROPEAN COMMITTEE FOR IRON AND STEEL STANDARDIZATION
COMITÉ EUROPÉEN DE NORMALISATION DU FER ET DE L'ACIER
EUROPÄISCHES KOMITEE FÜR EISEN-UND STAHLNORMUNG
EUROPEAN CERTIFIED REFERENCE MATERIAL (EURONORM – CRM)
CERTIFICATE OF CHEMICAL ANALYSIS
EURONORM – CRM No. 379-1 High alloy steel

Similar to EN 1.4563
LABORATORY MEANS (4 values) - Mass content in %

Line No	C	Si	Mn	P	S	Cr	Mo	Ni	B	Co	Cu	N	Sn	V
1	0.0112	-----	-----	0.0149	0.0005	26.590	3.2199	30.702	0.00173	0.0354	0.9645	0.0522	0.0015	0.0616
2	0.0114	0.3813	1.7898	0.0152	0.0005	26.624	3.2602	30.705	0.00177	0.0354	0.9688	0.0539	0.0017	0.0624
3	0.0115	0.3827	1.7927	0.0157	0.0005	26.657	3.2705	30.716	0.00180	-----	0.9732	0.0544	0.0017	0.0642
4	0.0118	0.3850	1.7940	0.0157	0.0006	26.703	3.2775	30.728	0.00182	0.0375	0.9757	0.0546	0.0018	0.0651
5	0.0119	0.3858	1.7945	0.0160	0.0006	26.710	3.2792	30.750	0.00183	0.0376	0.9761	0.0546	0.0019	0.0653
6	0.0119	0.3864	1.7948	0.0164	0.0006	26.713	3.2810	30.768	0.00183	0.0377	0.9778	0.0546	0.0021	0.0654
7	0.0120	0.3870	-----	0.0166	0.0007	26.733	3.2815	30.780	0.00183	0.0382	0.9787	0.0548	0.0022	0.0657
8	0.0121	0.3890	1.7968	0.0167	0.0007	26.753	3.2873	30.863	0.00185	0.0385	0.9791	0.0549	0.0023	0.0663
9	0.0121	0.3900	1.8005	0.0167	0.0007	26.765	3.2968	30.885	0.00188	0.0388	0.9805	0.0551	0.0023	0.0667
10	0.0122	0.3926	1.8030	0.0168	-----	26.765	3.3047	30.901	0.00188	0.0392	0.9823	0.0555	0.0023	0.0668
11	0.0123	0.3932	1.8030	0.0174	0.0007	26.780	3.3083	30.903	0.00192	0.0393	-----	0.0555	0.0024	0.0675
12	0.0124	0.3953	1.8036	0.0177	0.0009	26.819	3.3087	30.905	0.00203	0.0397	0.9907	-----	0.0024	0.0676
13	0.0128	0.3965	1.8098	0.0178	-----	26.837	3.3153	30.927	0.00208	0.0398	0.9940	0.0565	0.0024	0.0680
14	-----	0.3988	1.8151	0.0181	-----	26.844	3.3153	30.950	0.00210	0.0402	0.9984	0.0565	0.0025	0.0683
15	0.0135	0.4010	1.8188	-----	-----	26.844	3.3208	30.968	0.00218	0.0410	0.9995	0.0570	0.0069	0.0699
16	-----	0.4034	1.8212	-----	-----	26.888	3.3467	-----	0.0425	0.0410	1.0008	0.0412	0.0035	0.0008
17	-----	0.4056	1.8216	-----	-----	26.966	-----	-----	0.0425	0.0412	1.0035	-----	0.0008	0.0008
18	-----	0.4098	-----	-----	-----	26.975	-----	-----	0.0425	0.0412	1.0035	-----	0.0008	0.0008
19	-----	-----	-----	-----	-----	26.988	-----	-----	0.0425	0.0412	1.0035	-----	0.0008	0.0008
20	-----	-----	-----	-----	-----	-----	-----	-----	0.0425	0.0412	1.0035	-----	0.0008	0.0008
M_M	0.0121	0.3931	1.8039	0.0166	0.0006	26.787	3.2902	30.830	0.00190	0.0390	0.9840	0.0550	0.0021	0.0663
s_M	0.0006	0.0084	0.0109	0.0010	0.0002	0.114	0.0294	0.097	0.00014	0.0020	0.0122	0.0013	0.0004	0.0024
s_w	0.0003	0.0042	0.0078	0.0003	0.0001	0.046	0.0164	0.059	0.00012	0.0005	0.0041	0.0005	0.0002	0.0006

Line No	Ca	Sb	Al	As	Bi	Fe	O	Nb	Pb	Ti
1	0.0025	0.00048	0.0204	0.0018	0.000002	-----	0.0017	0.0012	0.000015	0.0008
2	0.0025	0.00053	0.0210	0.0026	0.000004	35.194	0.0018	0.0015	0.000021	0.0008
3	0.0027	0.00053	0.0222	0.0026	0.000005	35.417	0.0022	0.0017	-----	0.0010
4	0.0029	0.00054	0.0228	0.0030	0.000010	35.463	0.0023	0.0019	0.000025	0.0010
5	0.0029	0.00058	0.0229	0.0033	< 0.00002	35.557	0.0025	0.0020	0.000030	0.0011
6	0.0031	0.00058	0.0242	-----	≤ 0.00002	35.646	0.0028	0.0020	0.000100	0.0013
7	0.0032	0.00061	0.0250	-----	-----	35.705	0.0030	0.0039	-----	0.0022
8	0.0033	0.00062	0.0250	-----	-----	35.710	0.0034	0.0040	-----	0.0022
9	0.0033	0.00062	0.0267	-----	-----	35.765	0.0036	0.0041	0.00053	0.0026
10	0.0034	0.00064	0.0270	-----	-----	35.972	0.0041	-----	-----	-----
11	0.0036	-----	0.0273	-----	-----	-----	-----	-----	-----	-----
12	0.0037	-----	0.0277	-----	-----	-----	-----	-----	-----	-----
13	0.0042	-----	0.0281	-----	-----	-----	-----	-----	-----	-----
M_M	0.0033	0.00057								
s_M	0.0006	0.00006								
s_w	0.0002	0.00003								

M_M: Mean of the intralaboratory means
s_M: Standard deviation of the intralaboratory means
s_w: Intralaboratory standard deviation

The laboratory mean values have been examined statistically to eliminate outlying values.
Where a "-----" appears in the table it indicates that an outlying value has been eliminated by either the Cochran or Grubbs Test.

Values given in *italics* are for

Additional values from laboratories, for information (in µg/g):
Ag: 0.7, Ce: 0.1, Cs: 0.01, Ga: 17 and 29, Ir: 0.1, Mg: 6, Nd: 0.6, Os: 0.07, Pr: 0.2, Pt: 0.04, Rb: 4.3, Re: 2.4, Rh: 0.4, Ru: 1.4, Sm: 0.6, Ta: 0.04, W: 82, 92 and 100, Y: 0.1, Zr: 2.5 and 4.0

CERTIFIED VALUES - Mass content in %

	C	Si	Mn	P	S	Cr	Mo	Ni
M_M	0.0121	0.393	1.804	0.0166	0.0006	26.79	3.290	30.83
C (95 %)	0.0004	0.005	0.007	0.0006	0.0001	0.06	0.016	0.06
	B	Co	Cu	N	Sn	V	Ca	Sb
M_M	0.00190	0.0390	0.984	0.0550	0.0021	0.0663	0.0033	0.00057
C (95 %)	0.00008	0.0010	0.007	0.0007	0.0002	0.0013	0.0004	0.00004

The half-width confidence interval C(95%) = $\frac{t \times s_M}{\sqrt{n}}$ where "t" is the appropriate Student's t value and "n" is the number of acceptable mean values

For further information regarding the confidence interval for the certified value see ISO Guide 35:2006 sections 6.1 and 10.5.2

This reference material was prepared in accordance with the recommendations set out in
ISO Guides 30 – 35 and issued by:

Swerea | KIMAB

Drottning Kristinas väg 48, SE 114 28, Stockholm, Sweden

On behalf of: The Iron and Steel Nomenclature Co-ordinating Committee (COCOR) of the ECIISS,
after approval by all the participating laboratories and all the producing organisations.

(France-IRSID/CTIF; Germany-Iron and Steel CRM Working Group: Stahlinstitut VDEh,
BAM Bundesanstalt für Materialforschung und -prüfung & MPI für Eisenforschung;

UK-BAS Ltd; Nordic Countries-Nordic CRM Working Group)



SEPTEMBER 2010

EURONORM – CRM No. 379-1
METHODS USED

Element	Line number	Analytical methods
C	1.2.3.4.5.6.7.8.9.10.11.12.13 15	Combustion - infrared absorption Combustion - non-aqueous titration after absorption in organic solvent
Si	2.5.6.7.8.9.10.12.17 3.4.11.13.14.15.16.18	Gravimetry, dehydration with perchloric acid ICP-OES
Mn	2.8.9.17 3.5.6.10.11.13.15.16 4 12 14	Spectrophotometry, periodate oxidation ICP-OES PAA FAAS XRF
P	1.2.6.8.11.14 3.4.5.9.12.13 7.10	ICP-OES Spectrophotometry, phosphovanadomolybdate, extraction Spectrophotometry, molybdenum blue, extraction
S	1.2.3.4.5.6.7.8.9.11.12	Combustion - Infrared absorption
Cr	1.2.3.5.7.8.9.11.13.14.15.16 4.12 6.17.18.19 10	Titration with Fe (II), oxidation with persulphate Titration with Fe (II), oxidation with peroxide ICP-OES XRF
Mo	1 2.4.5.6.8.9.11.12.13.14 3.16 7 10	Spectrophotometry, thiocyanate in presence of Sn (II), without extraction ICP-OES Spectrophotometry, thiocyanate in presence of Sn (II), extraction XRF Gravimetry with 8-hydroxyquinoline, after anion-exchange
Ni	1.4.6.8.9.11.12.15 2 3.5.7.14 10 13	Gravimetry, dimethylglyoxime Titration with dichromate, separation with dimethylglyoxime ICP-OES Cyanometric titration Complexometric titration, visual end point
B	1.3.5.8.11 2.10 4.6.7.9.12.13.14.15	Spectrophotometry, curcumin ICP-MS ICP-OES
Co	1.4.6.7.8.9.11.12.15.17.18 2.5.10 13.14 16	ICP-OES FAAS NAA ETAAS
Cu	1.5.6.7.8.9.10.15.16.17 2.4.12.13.14 3	ICP-OES FAAS Spectrophotometry, diethyldithiocarbamate, extraction
N	1.2.3.4.5.6.7.8.9.10.11.14.15 12 13	Thermal conductivity, decomposition in graphite crucible Acidimetric titration after distillation, visual end point Spectrophotometry, indophenol blue, distillation
Sn	1.2.4.5.10.12 3.7.8.11.13.14 6 9	ICP-OES ETAAS FAAS, extraction with TOPO/KI/MIBK ICP-MS
V	1.3.4.6.8.9.11.12.13.14.15.16 2.5 7 10	ICP-OES FAAS NAA ETAAS
Ca	1.2.5.8.11.13 3.4.6.7.9.10.12.14	FAAS ICP-OES
Sb	1.2.3.9.10 4.6.7.8 5	ETAAS ICP-MS NAA

Abbreviations:

ETAAS	Electrothermal Atomic Absorption Spectrometry
FAAS	Flame Atomic Absorption Spectrometry
ICP-OES	Inductively Coupled Plasma – Optical Emission Spectrometry
ICP-MS	Inductively Coupled Plasma – Mass Spectrometry
NAA	Neutron Activation Analysis
PAA	Photon Activation Analysis
XRF	X-ray Fluorescence Spectrometry

<i>Element</i>	<i>Line number</i>	<i>Analytical methods</i>
<i>Al</i>	1.2.3.5.6.7.10.13 4.8.9.11.12	<i>ICP-OES</i> <i>FAAS</i>
<i>As</i>	1 2.5 3 4	<i>ICP-MS</i> <i>ETAAS</i> <i>FAAS</i> <i>NAA</i>
<i>Bi</i>	1.2.4.5.6 3	<i>ETAAS</i> <i>AAS, hydride</i>
<i>Fe</i>	2 3.4.5.7.9 6 8 10	<i>Titration with Cr (VI) after reduction with Sn (II)</i> <i>ICP-OES</i> <i>Titration with Mn (VII) after reduction with Ti (III)</i> <i>Titration with Mn (VII) after reduction with Sn (II)</i> <i>XRF</i>
<i>O</i>	1 2.3.4.5.6.7.8.9.10	<i>Thermal conductivity, reduction fusion under helium</i> <i>Infrared absorption measurement, reduction fusion under helium</i>
<i>Nb</i>	1.2.3.4.5.6.7.8.9.10	<i>ICP-OES</i>
<i>Pb</i>	1.2.4.5.6	<i>ETAAS</i>
<i>Ti</i>	1.2.3.4.5.6.7.9 8	<i>ICP-OES</i> <i>ETAAS</i>

Abbreviations:

ETAAS	Electrothermal Atomic Absorption Spectrometry
FAAS	Flame Atomic Absorption Spectrometry
ICP-OES	Inductively Coupled Plasma – Optical Emission Spectrometry
ICP-MS	Inductively Coupled Plasma – Mass Spectrometry
NAA	Neutron Activation Analysis
XRF	X-ray Fluorescence Spectrometry

DESCRIPTION OF THE SAMPLE

ECRM 379-1 is available in the form of milling chips in bottles containing 100 g. It is also available as 38 mm diameter discs 25 mm thick. The chips were passed through a 2000 µm aperture sieve and further sieving was carried out to exclude chips passing through a 250 µm aperture sieve.

INTENDED USE & STABILITY

The chip sample, ECRM 379-1(C), is intended for the verification of analytical methods, such as those used by the participating laboratories, for the calibration of analytical instruments in cases where the calibration with primary substances (pure stoichiometric metals or compounds) is not possible and for establishing values for secondary reference materials.

It will remain stable provided that the bottle remains sealed and is stored in a cool, dry atmosphere. When the bottle has been opened the lid should be secured immediately after use. If the contents should become discoloured (eg oxidised) due to atmospheric contamination they should be discarded.

The solid (disc) sample, ECRM 379-1(D), is intended for establishing and checking the calibration of instruments, such as Optical Emission Spectrometers and X-ray Spectrometers, for the analysis of samples of similar materials. The "as received" working surface of the sample should be finished before use in order to remove any protective coating. It will remain stable provided that it is not subjected to excessive heat (eg, during preparation of the working surface).

TRACEABILITY

The traceability of ECRM 379-1 has been established in accordance with principles of ISO Guides 30 – 35 and the International vocabulary of basic and general terms in metrology.

The characterisation of this material has been achieved by inter-laboratory study, each laboratory using the method of their choice, details of which are given above. These methods are either stoichiometric analytical techniques or methods which are calibrated against pure metals or stoichiometric compounds. Most methods used were either international or national standard methods or methods which are technically equivalent.

PARTICIPATING LABORATORIES

AB Sandvik Materials Technology, Sandviken (Sweden)
 Aciéries Aubert et Duval, Les Ancizes (France)
 AG der Dillinger Hüttenwerke, Dillingen (Germany)
 ALS Scandinavia AB, Luleå (Sweden)
 ArcelorMittal Imphy, Imphy (France)
 ArcelorMittal Isbergues Stainless, Isbergues (France)
 BAM Bundesanstalt für Materialforschung und -prüfung, Berlin (Germany)
 Bodycote Materials Testing, Teesside (UK)
 Böhler Edelstahl GmbH, Kapfenberg (Austria)
 CTIF, Charleville (France)
 CTIF, Sèvres (France)
 Kanthal AB, Hallstahammar (Sweden)
 Max-Planck-Institut für Eisenforschung, Düsseldorf (Germany)
 Outokumpu Stainless, Avesta (Sweden)
 Ovako Hofors AB, Hofors, (Sweden)
 Ridsdale & Co Ltd, Middlesbrough (UK)
 SCK-CEN, Belgian Nuclear Research Centre, Mol (Belgium)
 Sheffield Testing Laboratories, Sheffield (UK)
 SSAB, Oxelösund (Sweden)
 Swerea KIMAB AB, Stockholm (Sweden)
 Tata Steel Speciality, Rotherham (UK)
 ThyssenKrupp Nirosta, Krefeld (Germany)
 ThyssenKrupp VDM, Werdohl (Germany)
 Trinec Iron and Steel Works, Chemical and Mechanical Testing Laboratories, Třinec - Staré Město (Czech Republic)
 voestalpine Stahl GmbH, Linz (Austria)

EURONORM – CRM No. 379-1

FURTHER INFORMATION

For information regarding the preparation, certification and supply of these European Certified Reference Materials (EURONORM-CRMs) and the use of the statistical information given on this certificate, please refer either to the producer of this Certified Reference Material or to Technical Reports CEN/TR 10317:2009 and CEN/TR 10350:2009, both of which are available from the national standards body in your country. (In the UK this is the BSI, 389 Chiswick High Road, London W4 4AL).

Further information and advice on this or other Certified Reference Materials or Reference Materials produced by Swerea KIMAB AB, may be obtained from the address below.

Angaben über Herstellung, Zertifizierung und Bezugsmöglichkeiten dieser Europäischen Zertifizierten Referenzmaterialien (EURONORM-ZRM) sowie über die Anwendungen der in diesem Zertifikat enthaltenen statistischen Daten sind erhältlich beim Hersteller dieses zertifizierten Referenzmaterials, dessen Adresse auf diesem Zertifikat angegeben ist oder sie finden sich in den CEN-Reports CEN/TR 10317 und CEN/TR 10350, beide zu beziehen durch die nationalen Normenorganisationen (in Deutschland bei der Vertriebsstelle des DIN: Beuth-Verlag GmbH, Burggrafenstraße 4-10, 10787 Berlin).

Weitere Informationen und Hinweise zu diesem oder anderen durch Swerea KIMAB AB, hergestellten zertifizierten Referenzmaterialien oder Referenzmaterialien könnten unter der unten angegebenen Adresse erhalten werden.

Pour disposer d'informations sur la fabrication, la certification et la distribution des Matériaux de Référence Certifiés Européens (EURONORM-MRC) ainsi que sur l'utilisation des informations statistiques données sur ce certificat, se reporter soit au producteur de ce Matériau de Référence Certifié, soit aux Rapports Techniques CEN/TR 10317:2009 et CEN/TR 10350:2009. On peut se procurer ces deux documents auprès des organismes nationaux de normalisation. (Pour la France: AFNOR, 11 Avenue Francis de Pressensé, 93571 – St Denis la Plaine Cedex).

D'autres informations et avis au sujet de ce Matériau de Référence Certifié, ou de tout autre Matériau de Référence Certifié ou Matériau de Référence produit par le Swerea KIMAB AB, peuvent être demandés en contactant l'adresse figurant dans le bas de ce Certificat.

För information angående tillverkning, certifiering och anskaffning av dessa europeiska certifierade referensmaterial (EURONORM CRM) och för användning av statistisk information, som angivits i detta certifikat, refereras antingen till producenten av detta certifierade referensmaterial eller till Teknisk Rapport CEN/TR 10317:2009 och CEN/TR 10350:2009 som kan erhållas från den nationella standardiseringssorganisationen. (Sverige: SIS, S:t Paulsgatan 6, SE-118 80 Stockholm, Finland: SFS, PL. 116, FIN-002 41, Helsingfors, Danmark: DS, Kollegievej 6, DK-Charlottenlund 2920, Norge: NSF, Drammensveien, 145 A , Postboks 353 Skøyen, NO-0213 Oslo, Island: STRI, Holtagardar, IS-104 Reykjavik).

Ytterligare information och rådfrågan om detta eller andra Certifierade Referensmaterial/Referensmaterial, producerade av Swerea KIMAB AB, kan erhållas från angiven adress enligt nedan.

Rein Vainik
Research leader in Analytical Chemistry

Swerea KIMAB AB
Drottning Kristinas väg 48
SE 114 28 Stockholm, Sweden
E-mail: crm@swerea.se
Website: www.swreakimab.se